

Partial Denture Solutions

ROE Dental Laboratory offers a diverse selection of partial denture solutions. Having fabricated over 73,000 partial dentures in the past 9 years, our experienced technical team brings deep industry knowledge to address the unique needs and preferences of your patients.

Traditional Partial

Traditional removable partial dentures use metal frameworks with acrylic bases to restore partially edentulous patients. ROE offers two material options for the framework: SLM printed metal or titanium.

Traditional Partial Denture - CoCr	\$466	12 days in lab
Traditional Partial Denture - Titanium	\$514	13 days in lab
Frame Only - CoCr	\$205	4 days in lab
Frame Only - Titanium	\$250	5 days in lab



Thermoplastic Partial

Thermoplastic partial dentures are ideal for patients with metal or acrylic allergies or for those looking for a clinically-unbreakable, esthetic solution. Thermoplastic Hybrids include a metal (CoCr) substructure added to aid in stability.


Thermoplastic set & finish complete (includes reset)	\$308	8 days in lab
Thermoplastic Hybrid (includes reset) - CoCr	\$413	12 days in lab
Thermoplastic Unilateral or Flipper (includes reset)	\$201	4 days in lab



Acetal Partial

Acetal resin partials are very biocompatible, making it an excellent choice for patients with sensitivities to metal or acrylic.

Acetal Partial Denture	\$461	12 days in lab
Acetal Frame Only	\$199	4 days in lab



All-Acrylic Partial (Flipper)

Acrylic partial dentures are made to be supported by the mucosa and are often are fabricated for temporary use. Requires clasps for retention.

All-Acrylic Partial Denture (Flipper)	\$152	5 days in lab
---------------------------------------	-------	---------------



Order Today www.roedentallab.com/Removables